

After three years of continuous testing, modification, refinement, and improvement, DAG is proud to present it has taken great effort to achieve the highest quality requirements. We believe that our dermal hydrator is the best in class, and we are confident that it meets the highest standard of quality, durability and versatility.



La Racine ADH Full-face Drug Assisted Treatment

New Botulinum Toxin Full Face Rephrase Treatment

Medical aesthetic technology, breakthrough, 32G-9-pin needles, no leakage

FULL-FACE WRINKLE REMOVAL FACE LIFTING SWEAT ODOR REMOVAL

La Racine ADH latest advancement in minimally invasive injection instrument

Specifically suitable for: "Botulinum Toxin injection", "hair regeneration", "anti-aging for face" and other drug-assisted treatments, the operation can achieve three precise control: "depth", "dose", "no leakage"

Depth of Injection

Dose control

Inject without leaking

By using ADH instrument with Botulinum Toxin, it can achieve a significant full-face lifting effect by injecting the dermis layer. Usage: With precision consistency in the delivery of the drug, it has significant effect on wrinkle removal, face lifting, sweat removal and body odor removal.

Product name	La Racine Automated Dermal Hydrator
Certification	ISO 10993 Biocompatibility test IEC 60601-1 Medical Electronics Safety / EMC IEC 60601-1-6 Usability test IEC 62304 Software validation Animal experiment
Operation mode	Three modes: discontinuous, continuous, fast
Needle type	Single needle, 9-pin needles
Weight	Lightweight (only 642g)
Wire	Medical grade TPE material, life expectancy increased by three times compared to normal wire, and will not decompose toxic substances
Filter	Patented filter tanks to achieve complete container cleaning and to prevent from the risk of blood contamination
Operation noise	Low noise (52-54dB), conforms to human body standards
Data storage	Wireless network storage of medical information for data management
Injection depth accuracy	±0.1 mm (The highest standard in the industry)
Injection dose accuracy	±2 ul (The highest standard in the industry, allows 120 rounds of injection for effective treatment)
Power supply	Meets IEC 60601-1 medical grade standards, safe and stable
Price	\$11,500 U.S.



DAG Regenerative Engineering Pte Ltd

Nordcom One @ Gambas, 3 Gambas Crescent #06-15 Singapore 757088
TEL : +65 8748 8183
URL : <http://www.dag-regen.com>
info@dag-regen.com

La Racine™

ADH

AUTOMATED DERMAL HYDRATOR

Third-Generation, Innovative Automated Dermal Hydrator

DEPTH CONTROL PRECISE DOSE NO LEAKAGE MINIMAL INVASION SWISS MADE MOTOR



La Racine ADH is the latest generation of minimally invasive cosmetic medical instrument for auxiliary injection. Using the patented designed gun to control the depth of injection with automatic negative-pressured suction technology, the agent is injected into the skin through multiple needles. Injection frequency, formula and the dosage can be adjusted which lead to a more precise aesthetic application of medicine, leading to optimal treatment outcome with high patient satisfaction.